

Anti-Rip3 Rabbit Polyclonal Antibody For Ihc-P, If/Icc, Elisa - Q9QzL0

Item Number: CM0020566



Introduction

Rabbit polyclonal antibody against mouse RIP3 (Ripk3), validated for IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Suitable for necroptosis, apoptosis, and cell death research.

[Learn More](#)

Parameter	Value
Product Name	RIP3 Rabbit pAb
Remarks/Alias	Rip3; 2610528K09Rik; RIP3
Species	Mouse
GeneID (Human)	11035
GeneID	56532
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 320-420 of mouse RIP3 (NP_064339.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9QZL0
Protein Weight	53kDa
Shipping	Ice bag