


Anti-Ripk1/Rip Rabbit Polyclonal Antibody For Wb, Elisa - Q13546

Item Number: CM0025634



Introduction

Validated rabbit polyclonal antibody targeting RIPK1/RIP (RIP1) for Western blot (WB) and ELISA. Detects human and mouse RIP1, a serine-threonine kinase regulating TNF-induced apoptosis, necroptosis, and inflammatory pathways. Suitable for cell death signaling research.

[Learn More](#)

Parameter	Value
Product Name	RIPK1/RIP Rabbit pAb
Remarks/Alias	RIP; RIP1; AIEFL; IMD57; RIP-1
Species	Human
GeneID	8737
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 465-671 of human RIPK1/RIP (NP_003795.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13546
Protein Weight	75kDa
Shipping	Ice bag