

Anti-Ripk3 (Rip3) Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Y572

Item Number: CM0005985



Introduction

Monoclonal antibody targeting human RIPK3 (RIP3), suitable for Western blot, IHC-P, and ELISA; cross-reacts with mouse and rat. Ideal for studying necroptosis and apoptosis pathways in cancer and viral infection research.

[Learn More](#)

Parameter	Value
Product Name	Ripk3 Mouse mAb
Remarks/Alias	RIP3; RIPK3
Species	Human
GeneID (human)	11035
GeneID	11035
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 101-200 of human RIP3 (NP_006862.2).
Source	Mouse
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y572
Protein Weight	57kDa
Shipping	Ice bag