


Anti-RLN2 Polyclonal Antibody For Wb, Ihc-P, Elisa - P04090

Item Number: CM0025092



Introduction

Rabbit polyclonal antibody targeting human prorelaxin H2 (RLN2), validated for WB, IHC-P, and ELISA. Reacts with human and rat. Ideal for studying pregnancy-related connective tissue remodeling and reproductive biology.

[Learn More](#)

Parameter	Value
Product Name	RLN2 Rabbit pAb
Remarks/Alias	H2; RLXH2; H2-RLX; bA12D24.1.1; bA12D24.1.2; RLN2
Species	Human
GeneID (Human)	6019
GeneID	6019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-185 of human RLN2 (NP_604390.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P04090
Protein Weight	21kDa
Shipping	Ice bag