

Anti-Rnase13 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q5Gan3

Item Number: CM0015506

Introduction

Rabbit polyclonal antibody against human RNASE13, validated for WB, IHC-P, ELISA. Detects endogenous levels in human and mouse. Ideal for studying secreted inactive ribonuclease-like proteins.

[Learn More](#)

Parameter	Value
Product Name	RNASE13 Rabbit pAb
Remarks/Alias	RAL1; HEL-S-86p; RNASE13
Species	Human
GeneID (Human)	440163
GeneID	440163
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human RNASE13 (NP_001012264.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5GAN3
Protein Weight	18kDa
Shipping	Ice bag