

Anti-Rnf11 Polyclonal Antibody For Wb, Elisa - Q9Y3C5

Item Number: CM0029016

Introduction

Rabbit polyclonal antibody targeting human RNF11, validated for Western blot (WB) and ELISA, with cross-reactivity in mouse and rat. Ideal for studying ubiquitin-mediated inflammatory signaling and protein degradation pathways.

[Learn More](#)

Parameter	Value
Product Name	RNF11 Rabbit pAb
Remarks/Alias	CGI-123; SID1669; RNF11
Species	Human
GeneID (Human)	26994
GeneID	26994
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-154 of human RNF11 (NP_055187.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9Y3C5
Protein Weight	17kDa
Shipping	Ice bag