

Anti-Rnf166 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q96A37

Item Number: CM0024785

Introduction

Rabbit polyclonal antibody against human RNF166, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Ideal for studying E3 ubiquitin ligase function in innate immunity and autophagy.

[Learn More](#)

Parameter	Value
Product Name	RNF166 Rabbit pAb
Remarks/Alias	RNF166
Species	Human
GeneID (Human)	115992
GeneID	115992
Immunogen	Synthetic peptide, corresponding to amino acids 102-201 of human RNF166 (NP_849163.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q96A37
Protein Weight	26kDa
Shipping	Ice bag