


Anti-Rpb3/Polr2C Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P19387

Item Number: CM0020088



Introduction

Rabbit polyclonal antibody against DNA-directed RNA polymerase II subunit RPB3 (POLR2C). Validated for WB, IP, ELISA in human and mouse samples. Suitable for transcription and RNA polymerase II research.

[Learn More](#)

Parameter	Value
Product Name	RPB3/POLR2C Rabbit pAb
Remarks/Alias	RPB3; RPB31; hRPB33; hSRPB3; RPB3/POLR2C
Species	Human
Gene ID (Human)	5432
Gene ID	5432
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-275 of human RPB3/POLR2C (NP_116558.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P19387
Protein Weight	31kDa
Shipping	Ice bag