


Anti-Rpl14 Polyclonal Antibody For Wb, If/Icc, Elisa - P50914

Item Number: CM0027687



Introduction

Rabbit anti-RPL14 polyclonal antibody for Western blot, IF/ICC, and ELISA. Recognizes human, mouse, rat. 23kDa ribosomal protein eL14, cytoplasm. Recombinant immunogen. Ideal for protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPL14 Rabbit pAb
Remarks/Alias	L14; RL14; eL14; hRL14; CTG-B33; CAG-ISL-7; RPL14
Species	Human
GeneID (Human)	9045
GeneID	9045
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-215 of human RPL14 (NP_003964.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50914
Protein Weight	23kDa
Shipping	Ice bag