


Anti-Rpl38 Polyclonal Antibody For Wb, Elisa - P63173

Item Number: CM0017150



Introduction

High-quality rabbit polyclonal antibody specific to human ribosomal protein L38 (RPL38). Validated for WB and ELISA, detects 8kDa band. Essential for studies on ribosome structure, translation, and protein synthesis.

[Learn More](#)

Parameter	Value
Product Name	RPL38 Rabbit pAb
Remarks/Alias	L38; eL38; RPL38
Species	Human
Gene ID	6169
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human RPL38 (NP_000990.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P63173
Protein Weight	8kDa
Shipping	Ice bag