


Anti-Rplp0 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P05388

Item Number: CM0029639



Introduction

High-quality rabbit polyclonal antibody against human RPLP0 (UniProt P05388). Validated for WB, IF/ICC, ELISA, with cross-reactivity to mouse. Ideal for ribosomal protein expression studies and loading control applications.

[Learn More](#)

Parameter	Value
Product Name	RPLP0 Rabbit pAb
Remarks/Alias	P0; LP0; L10E; RPP0; uL10; PRLP0; RPLP0
Species	Human
GeneID (human)	6175
GeneID	6175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-317 of human RPLP0 (NP_000993.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05388
Protein Weight	34kDa
Shipping	Ice bag