

Anti-Rpn1 Polyclonal Antibody For Wb, Elisa - P04843

Item Number: CM0027459



Introduction

Rabbit polyclonal antibody against human RPN1 (ribophorin I), a subunit of the oligosaccharyl transferase complex essential for protein N-glycosylation. Suitable for WB and ELISA, cross-reacts with human and rat.

[Learn More](#)

Parameter	Value
Product Name	RPN1 Rabbit pAb
Remarks/Alias	OST1; RBPH1; RPN1
Species	Human
Gene ID (Human)	6184
Gene ID	6184
Immunogen	Recombinant protein / recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human RPN1 (NP_002941.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P04843
Protein Weight	69kDa
Shipping	Ice bag