

Anti-Rpn1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P04843

Item Number: CM0015071



Introduction

Rabbit polyclonal antibody against human RPN1 (Ribophorin I), a 69 kDa OST complex subunit. Validated for WB, IHC-P, ELISA; cross-reacts with mouse and rat. Suitable for N-glycosylation and ER studies.

[Learn More](#)

Parameter	Value
Product Name	RPN1 Rabbit pAb
Remarks/Alias	OST1; RBPH1; RPN1
Species	Human
Gene ID (Human)	6184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-280 of human RPN1 (NP_002941.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P04843
Protein Weight	69kDa
Shipping	Ice bag