

Anti-Rpp25L Rabbit Polyclonal Antibody For Wb, Elisa - Q8N5L8

Item Number: CM0025320

Introduction

Rabbit polyclonal antibody against human RPP25L (C9orf23), a putative RNase P/MRP subunit. Validated for Western Blot and ELISA, cross-reacts with mouse. Suitable for nucleus/RNA processing research.

[Learn More](#)

Parameter	Value
Product Name	RPP25L Rabbit pAb
Remarks/Alias	C9orf23; bA296L22.5; RPP25L
Species	Human
GeneID (human)	138716
GeneID	138716
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-163 of human RPP25L (NP_680544.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q8N5L8
Protein Weight	18kDa
Shipping	Ice bag