

Anti-Rps18 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62269

Item Number: CM0012907



Introduction

Rabbit polyclonal antibody against human RPS18 (40S ribosomal protein S18), validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Suitable for ribosome research and protein synthesis studies.

[Learn More](#)

Parameter	Value
Product Name	RPS18 Rabbit pAb
Remarks/Alias	KE3, S18, HKE3, KE-3, uS13, D6S218E, RPS18
Species	Human
GeneID (Human)	6222
GeneID	6222
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human RPS18 (NP_072045.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62269
Protein Weight	18kDa
Shipping	Ice bag