


Anti-Rps6Ka3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P51812

Item Number: CM0025437



Introduction

Rabbit polyclonal antibody against human RPS6KA3 (RSK2, p90RSK2) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Detects ~84 kDa band. For ERK/MAPK signaling, translation, and cell survival studies.

[Learn More](#)

Parameter	Value
Product Name	RPS6KA3 Rabbit pAb
Remarks/Alias	CLS; RSK; HU-3; RSK2; MRX19; ISPK-1; XLID19; p90-RSK2; pp90RSK2; MAPKAPK1B; S6K-alpha3; RPS6KA3
Species	Human
Gene ID	6197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 68-327 of human RPS6KA3 (NP_004577.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P51812
Protein Weight	84kDa
Shipping	Ice bag