

Anti-Rps7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62081

Item Number: CM0029585



Introduction

High-quality rabbit polyclonal antibody against human ribosomal protein S7 (RPS7), validated for WB, IHC-P, IF/ICC, and ELISA applications; reacts with human, mouse, and rat samples. Ideal for ribosomal biogenesis and Diamond-Blackfan anemia research.

[Learn More](#)

Parameter	Value
Product Name	RPS7 Rabbit pAb
Remarks/Alias	S7; eS7; DBA8; RPS7
Species	Human
Gene ID (Human)	6201
Gene ID	6201
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-194 of human RPS7 (NP_001002.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62081
Protein Weight	22kDa
Shipping	Ice bag