

Anti-Rps8 Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P62241

Item Number: CM0003551



Introduction

High-quality rabbit monoclonal antibody against human RPS8 (40S ribosomal protein S8, SWISS: P62241). Validated for WB, IF/ICC, IP, ELISA with cross-reactivity to mouse and rat. Ideal for ribosomal research.

[Learn More](#)

Parameter	Value
Product Name	RPS8 Rabbit mAb
Remarks/Alias	S8; eS8; RPS8
Species	Human
Gene ID	6202
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-208 of human RPS8 (NP_001003.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62241
Protein Weight	24 kDa
Shipping	Ice bag