


Anti-Rragd Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Nql2

Item Number: CM0030656



Introduction

Rabbit polyclonal antibody against human Ras-related GTP-binding protein D (RRAGD) validated for WB, IHC-P, ELISA; cross-reactive with mouse and rat; detects 46kDa; key regulator of mTORC1 signaling in response to amino acids.

[Learn More](#)

Parameter	Value
Product Name	RRAGD Rabbit pAb
Remarks/Alias	RAGD; HOMG7; bA11D8.2.1; RRAGD
Species	Human
Gene ID	58528
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 181-400 of human RRAGD (NP_067067.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NQL2
Protein Weight	46kDa
Shipping	Ice bag