

Anti-Ruvbl1/Tip49A/Pontin Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa, Chip - Q9Y265

Item Number: CM0004130



Introduction

High-quality recombinant rabbit monoclonal antibody targeting RUVBL1/TIP49A/PONTIN (Q9Y265). Validated for Western blot, IHC-P, ELISA, and ChIP; cross-reacts with human, mouse, and rat. Suitable for studying chromatin remodeling, DNA repair, and oncogenic signaling.

[Learn More](#)

Parameter	Value
Product Name	RUVBL1/TIP49A/PONTIN Rabbit mAb
Remarks/Alias	RVB1; TIH1; ECP54; TIP49; ECP-54; INO80H; NMP238; PONTIN; TIP49A; NMP 238; Pontin52; RUVBL1/TIP49A/PONTIN
Species	Human
Gene ID (Human)	8607
Gene ID	8607
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 350-456 of human RUVBL1/TIP49A/PONTIN (Q9Y265)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA, ChIP
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y265
Protein Weight	50kDa
Shipping	Ice bag