

Anti-S100A10 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P60903

Item Number: CM0018470



Introduction

S100A10 Rabbit Polyclonal Antibody validated for WB, IHC-P, IF/ICC, ELISA. Detects human, mouse, rat S100A10. Suitable for studies of ANXA2 interaction and protein phosphorylation regulation.

[Learn More](#)

Parameter	Value
Product Name	S100A10 Rabbit pAb
Remarks/Alias	42C; P11; p10; GP11; ANX2L; CAL1L; CLP11; Ca[1]; ANX2LG; S100A10
Species	Human
GeneID (Human)	6281
GeneID	6281
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-97 of human S100A10 (NP_002957.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P60903
Protein Weight	11kDa
Shipping	Ice bag