


Anti-Sars-Cov-2 N Protein Rabbit Polyclonal Antibody For Wb, Elisa - P59595

Item Number: CM0020939



Introduction

Rabbit polyclonal antibody targeting SARS-CoV-2 nucleoprotein (N protein). Validated for WB and ELISA. Detects a ~46 kDa protein. Suitable for viral assembly, host interaction, and diagnostic research.

[Learn More](#)

Parameter	Specification
Product Name	SARS-CoV-2 N Protein Rabbit pAb
Gene ID	1489678
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of coronavirus N Protein (NP_828858.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P59595
Protein Weight	46kDa
Shipping	Ice bag