


Anti-Sars-Cov-2 Nsp15 Rabbit Polyclonal Antibody For Wb, Elisa - P0Dtd1

Item Number: CM0024312



Introduction

Rabbit polyclonal antibody targeting SARS-CoV-2 NSP15 endoribonuclease, suitable for Western blot and ELISA. Detects viral NSP15 involved in RNA replication and immune evasion. Ideal for COVID-19 research.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV-2 NSP15 Rabbit pAb
Gene ID	43740578
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 6600-6700 of coronavirus NSP15 (YP_009724389.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV-2
SWISS	P0DTD1
Protein Weight	794 kDa
Shipping	Ice bag