

Anti-Sars-Cov Spike Polyclonal Antibody For Wb & Elisa - P59594

Item Number: CM0020972



Introduction

Rabbit polyclonal antibody targeting SARS-CoV Spike glycoprotein (E2). Validated for Western blot and ELISA. Cross-reacts with SARS-CoV. Ideal for studying viral entry, fusion, and host immune evasion mechanisms.

[Learn More](#)

Core Product Specifications and Parameters	
Product Name	SARS-CoV Spike Rabbit pAb
Remarks/Alias	E2; SARS-CoV Spike
Gene ID	1489668
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of coronavirus Spike (NP_828851.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV
SWISS	P59594
Protein Weight	139kDa
Shipping	Ice bag