

Anti-Sars-Cov Spike Rbd Polyclonal Antibody For Wb, Elisa - P59594

Item Number: CM0021814



Introduction

Rabbit polyclonal antibody specific to the receptor-binding domain (RBD) of SARS-CoV Spike glycoprotein. Validated for WB and ELISA, it serves as a key reagent for coronavirus research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	SARS-CoV Spike RBD Rabbit pAb
Remarks/Alias	E2; SARS-CoV Spike RBD
Gene ID	1489668
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of coronavirus Spike RBD (NP_828851.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	SARS-CoV
SWISS	P59594
Protein Weight	139kDa
Shipping	Ice bag