

Anti-Scamp2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O15127

Item Number: CM0019534

Introduction

High-quality anti-SCAMP2 rabbit polyclonal antibody validated for WB, IF/ICC, ELISA. Recognizes human, mouse, rat SCAMP2 (37 kDa). Ideal for studying secretory carrier-associated membrane protein 2 in post-Golgi recycling.

[Learn More](#)

Parameter	Value
Product Name	SCAMP2 Rabbit pAb
Remarks/Alias	SCAMP2
Species	Human
Gene ID	10066 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-152 of human SCAMP2 (NP_005688.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15127
Protein Weight	37kDa
Shipping	Ice bag