


Anti-Scn2A Polyclonal Antibody For Wb, Elisa - Q99250

Item Number: CM0018366



Introduction

Rabbit polyclonal antibody targeting human SCN2A (Nav1.2 sodium channel alpha subunit). Suitable for WB and ELISA; 228 kDa. Ideal for neuroscience and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	SCN2A Rabbit pAb
Remarks/Alias	EA9; HBA; NAC2; BFIC3; BFIS3; BFNIS; DEE11; HBSCI; EIEE11; HBSCII; Nav1.2; SCN2A1; SCN2A2; Na(v)1.2; SCN2A
Species	Human
GeneID (Human)	6326
GeneID	6326
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 275-401 of human SCN2A (NP_001035232.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q99250
Protein Weight	228kDa
Shipping	Ice bag