

Anti-Sdf4 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Brk5

Item Number: CM0017692

Introduction

Anti-SDF4 rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Detects human and mouse SDF4, a 45 kDa calcium-binding protein involved in calcium-dependent ER activities. Ideal for metabolic and Golgi research.

[Learn More](#)

Parameter	Value
Product Name	SDF4 Rabbit pAb
Remarks/Alias	Cab45; SDF-4; SDF4
Species	Human
GeneID (Human)	51150
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 37-200 of human SDF4 (NP_057260.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BRK5
Protein Weight	42kDa
Shipping	Ice bag