

Anti-Sdhb Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P21912

Item Number: CM0014449



Introduction

Rabbit polyclonal antibody against human SDHB (succinate dehydrogenase iron-sulfur subunit). Validated for Western blot (WB), immunofluorescence/immunocytochemistry (IF/ICC), and ELISA. Cross-reacts with mouse and rat SDHB. Ideal for mitochondrial complex II research, oxidative phosphorylation studies, and metabolic disease investigations.

[Learn More](#)

Parameter	Value
Product Name	SDHB Rabbit pAb
Remarks/Alias	IP; SDH; CWS2; PGL4; SDH1; SDH2; SDHIP; MC2DN4; SDHB
Species	Human
Gene ID (Human)	6390
Gene ID	6390
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-280 of human SDHB (NP_002991.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21912
Protein Weight	32kDa
Shipping	Ice bag