

Anti-Sec11A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P67812

Item Number: CM0016755



Introduction

High-quality rabbit polyclonal antibody against human SEC11A, validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Ideal for studying endoplasmic reticulum signal peptide cleavage.

[Learn More](#)

Parameter	Value
Product Name	SEC11A Rabbit pAb
Remarks/Alias	SPC18, SPCS4A, SEC11L1, sid2895, 1810012E07Rik, SEC11A
Species	Human
GeneID (Human)	23478
GeneID	23478
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 37-145 of human SEC11A (NP_055115.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67812
Protein Weight	21kDa
Shipping	Ice bag