

Anti-Sec13 Rabbit Polyclonal Antibody For Wb, Ihc, If/Icc, Elisa - P55735

Item Number: CM0014241



Introduction

High-quality rabbit polyclonal antibody against human SEC13, validated for WB, IHC, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Ideal for studying NPC, COPII, and mTORC1 signaling.

[Learn More](#)

Parameter	Specification
Product Name	SEC13 Rabbit pAb
Remarks/Alias	SEC13R; npp-20; SEC13L1; D3S1231E; SEC13
Species	Human
GeneID (Human)	6396
GeneID	6396
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-322 of human SEC13 (NP_899195.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P55735
Protein Weight	36kDa
Shipping	Ice bag