


Anti-Sel1L Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Ubv2

Item Number: CM0005833



Introduction

Rabbit monoclonal antibody against human SEL1L, validated for WB, IHC-P, and ELISA. Cross-reacts with mouse and rat. Targets ERAD and Notch signaling, ideal for pancreatic research.

[Learn More](#)

Parameter	Value
Product Name	SEL1L Rabbit mAb
Remarks/Alias	Hrd3; SEL1L1; PRO1063; SEL1-LIKE
Species	Human
Gene ID	6400
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 251-350 of human SEL1L (NP_005056.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UBV2
Protein Weight	89kDa
Shipping	Ice bag