

Anti-Septin 8 Rabbit Monoclonal Antibody For Wb And Elisa - Q92599

Item Number: CM0003593



Introduction

High-quality Septin 8 rabbit monoclonal antibody validated for WB and ELISA. Reacts with human, mouse, and rat Septin-8, a cytoskeletal GTPase involved in platelet secretion and BACE1 trafficking. Ideal for Alzheimer's and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	Septin 8 Rabbit mAb
Remarks/Alias	SEP2; SEPT8; Septin 8
Species	Human
Gene ID (Human)	23176
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human Septin 8 (NP_001287728.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92599
Protein Weight	56kDa
Shipping	Ice bag