

Anti-Septin 9 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uhd8

Item Number: CM0028823



Introduction

Highly specific rabbit polyclonal antibody targeting human Septin-9. Validated for WB, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat. Suitable for cytoskeletal dynamics, cytokinesis, and microbial pathogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Septin 9 Rabbit pAb
Remarks/Alias	MSF; MSF1; NAPB; SEPT9; SINT1; PNUTL4; SeptD1; AF17q25; Septin 9
Species	Human
Gene ID (Human)	10801
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human Septin 9 (NP_001106964.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UHD8
Protein Weight	65kDa
Shipping	Ice bag