

Anti-Serinc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, And Elisa - Q9Nrx5

Item Number: CM0021605

Introduction

Rabbit polyclonal antibody targeting human SERINC1 (serine incorporator 1), validated for Western blot, IHC-P, and ELISA applications, with cross-reactivity to mouse and rat. Suitable for investigating phosphatidylserine and sphingolipid biosynthesis.

[Learn More](#)

Parameter	Value
Product Name	SERINC1 Rabbit pAb
Remarks/Alias	TDE2; TMS2; TDE1L; TMS-2; SERINC1
Species	Human
GeneID (human)	57515
GeneID	57515
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 350-453 of human SERINC1 (NP_065806.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRX5
Protein Weight	50kDa
Shipping	Ice bag