


Anti-Serpine2 Polyclonal Antibody For Wb, Elisa - P07093

Item Number: CM0016119



Introduction

Rabbit polyclonal antibody targeting SERPINE2 (Glia-derived nexin/Protease nexin 1), validated for WB and ELISA with cross-reactivity to mouse and rat. Swiss P07093.

[Learn More](#)

Parameter	Value
Product Name	SERPINE2 Rabbit pAb
Remarks/Alias	GDN; PI7; PN1; PNI; PI-7; PN-1; GDNPF; SERPINE2
Species	Human
GeneID (Human)	5270
GeneID	5270
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P07093
Protein Weight	44kDa
Shipping	Ice bag