


Anti-Ser/Thr-Pp2A Activator (Ptpa) Polyclonal Antibody For Wb, If/Icc, Elisa - Q15257

Item Number: CM0022584



Introduction

High-quality rabbit polyclonal antibody against human PTPA/PPP2R4, validated for Western blot, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Ideal for investigating PP2A regulation, apoptosis, and neurodegenerative disease mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Ser/thr-PP2A activator (PTPA) Rabbit pAb
Remarks/Alias	PP2A; PR53; PPP2R4; Ser/thr-PP2A activator (PTPA)
Species	Human
GeneID (human)	5524
GeneID	5524
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-323 of human Ser/thr-PP2A activator (PTPA) (NP_821067.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15257
Protein Weight	41kDa
Shipping	Ice bag