

Anti-Setd5 Rabbit Polyclonal Antibody For If/Icc, Elisa - Q9C0A6

Item Number: CM0014484

Introduction

SETD5 rabbit polyclonal antibody validated for IF/ICC and ELISA. Detects human SETD5, a histone methyltransferase involved in brain development and chromatin regulation. Ideal for neurodevelopmental and epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	SETD5 Rabbit pAb
Remarks/Alias	MRD23; SETD5A; SETD5
Species	Human
GeneID (Human)	55209
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-740 of human SETD5 (NP_001073986.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9C0A6
Protein Weight	158kDa
Shipping	Ice bag