


Anti-Sf3A3 Rabbit Monoclonal Antibody For Wb, Elisa - Q12874

Item Number: CM0003944



Introduction

Anti-SF3A3 rabbit monoclonal antibody for WB and ELISA, recognizing human and mouse SF3A3 (59 kDa). Essential tool for spliceosome research and pre-mRNA splicing studies.

[Learn More](#)

Parameter	Value
Product Name	SF3A3 Rabbit mAb
Remarks/Alias	PRP9; PRPF9; SAP61; SF3a60; SF3A3
Species	Human
GeneID (Human)	10946
GeneID	10946
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SF3A3 (Q12874)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q12874
Protein Weight	59kDa
Shipping	Ice bag