


Anti-Sf3B4 Rabbit Polyclonal Antibody For Wb And Elisa - Q15427

Item Number: CM0023810



Introduction

Rabbit polyclonal anti-SF3B4 antibody validated for WB and ELISA with human, mouse, rat samples. Targets splicing factor 3B subunit 4 for spliceosome research - Q15427.

[Learn More](#)

Parameter	Value
Product Name	SF3B4 Rabbit pAb
Remarks/Alias	AFD1; Hsh49; SAP49; SF3b49; SF3B4
Species	Human
GeneID (Human)	10262
GeneID	10262
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human SF3B4 (NP_005841.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15427
Protein Weight	44kDa
Shipping	Ice bag