


Anti-Sftpc Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P11686

Item Number: CM0007062



Introduction

Rabbit monoclonal antibody against human surfactant protein C (SFTPC) validated for WB, IF-P, and ELISA. Cross-reacts with mouse and rat. Suitable for respiratory disease research and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SFTPC Rabbit mAb
Remarks/Alias	SP5; SP-C; PSP-C; SFTP2; SMDP2; BRICD6; SFTPC
Species	Human
Gene ID (Human)	6440
Immunogen	Recombinant protein of human surfactant protein c/SP-C.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P11686
Protein Weight	21kDa
Shipping	Ice bag