


Anti-Sgol1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q5Fbb7

Item Number: CM0022462



Introduction

SGOL1 Rabbit pAb (polyclonal, Rabbit, WB/IF/ICC/ELISA, Human/Mouse/Rat) targets Shugoshin 1 (Q5FBB7), a key mitotic cohesion regulator. For IVD research and assay development by CamelBio.

[Learn More](#)

Parameter	Value
Product Name	SGOL1 Rabbit pAb
Remarks/Alias	SGO; CAID; SGOL1; NY-BR-85
Species	Human
GeneID (Human)	151648
GeneID	151648
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 10-100 of human SGOL1 (NP_612493.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5FBB7
Protein Weight	64kDa
Shipping	Ice bag