

Anti-Shox Rabbit Polyclonal Antibody For Wb, Elisa - 015266

Item Number: CM0022434

Introduction

Rabbit polyclonal antibody against human SHOX (short stature homeobox protein). Suitable for Western blot and ELISA applications, cross-reacts with mouse and rat. Designed for growth and skeletal development studies.

[Learn More](#)

Parameter	Value
Product Name	SHOX Rabbit pAb
Remarks/Alias	SS; GCFX; PHOG; SHOXY; SHOX
Species	Human
Gene ID (Human)	6473
Gene ID	6473
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-110 of human SHOX (NP_000442.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15266
Protein Weight	32kDa
Shipping	Ice bag