


Anti-Shp2 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q06124

Item Number: CM0002934



Introduction

High-quality rabbit monoclonal antibody against human SHP2 (PTPN11), validated for WB, IF/ICC, and ELISA. Recognizes human, mouse, and rat SHP2, ideal for signal transduction research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	SHP2 Rabbit mAb
Remarks/Alias	CFC; NS1; JMML; SHP2; BPTP3; PTP2C; METCDS; PTP-1D; SH-PTP2; SH-PTP3
Species	Human
GeneID (Human)	5781
GeneID	5781
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SHP2 (Q06124).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q06124
Protein Weight	68kDa
Shipping	Ice bag