


Anti-Sigmar1 Rabbit Polyclonal Antibody For Wb And Elisa - Q99720

Item Number: CM0017391



Introduction

Rabbit polyclonal antibody against human SIGMAR1, validated for WB and ELISA with cross-reactivity to human and rat. Suitable for research on neurodegenerative diseases, calcium signaling, and lipid transport.

[Learn More](#)

Parameter	Value
Product Name	SIGMAR1 Rabbit pAb
Remarks/Alias	SRBP; ALS16; DSMA2; OPRS1; SR-BP; SIG-1R; SR-BP1; sigma1R; hSigmaR1; SIGMAR1
Species	Human
Gene ID (Human)	10280
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-80 of human SIGMAR1 (NP_005857.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q99720
Protein Weight	25kDa
Shipping	Ice bag