

Anti-Sim2 Rabbit Polyclonal Antibody For Wb, Elisa - Q14190

Item Number: CM0012919

Introduction

Validated for WB and ELISA in human, mouse, and rat. SIM2 is a bHLH transcription factor critical for CNS development. Ideal for neuroscience and developmental biology studies.

[Learn More](#)

Parameter	Value
Product Name	SIM2 Rabbit pAb
Remarks/Alias	SIM; bHLHe15; HMC13F06; HMC29C01; SIM2
Species	Human
Gene ID (Human)	6493
Gene ID	6493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 342-525 of human SIM2 (NP_005060.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14190
Protein Weight	73kDa
Shipping	Ice bag