

Anti-Sirt5 Rabbit Monoclonal Antibody For Wb And Elisa - Q9Nxa8

Item Number: CM0007663



Introduction

High-quality rabbit monoclonal antibody against human SIRT5 (SIR2L5), validated for WB and ELISA. Recognizes human, mouse, and rat SIRT5; ideal for studying deacylase activity and metabolic regulation.

[Learn More](#)

Parameter	Value
Product Name	SIRT5 Rabbit mAb
Remarks/Alias	SIR2L5; SIRT5
Species	Human
GeneID (Human)	23408
GeneID	23408
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 37-310 of human SIRT5 (NP_036373.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NXA8
Protein Weight	34kDa
Shipping	Ice bag