


Anti-Slc18A3 Rabbit Polyclonal Antibody For If-P, Ihc-P, Elisa - Q16572

Item Number: CM0020130



Introduction

Rabbit polyclonal antibody against SLC18A3 (VACHT), validated for IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for cholinergic neurotransmission and synaptic vesicle research.

[Learn More](#)

Parameter	Value
Product Name	SLC18A3 Rabbit pAb
Remarks/Alias	CMS21; VACHT; SLC18A3
Species	Human
GeneID (human)	6572
GeneID	6572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 473-532 of human SLC18A3 (NP_003046.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q16572
Protein Weight	57kDa
Shipping	Ice bag