


Anti-Slc1A1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P43005

Item Number: CM0030452



Introduction

Rabbit polyclonal anti-SLC1A1/EAAT3 antibody validated for Western blot, immunofluorescence (IF-P), and ELISA. Detects human, mouse, and rat excitatory amino acid transporter 3. For neuroscience, renal, and oxidative stress research.

[Learn More](#)

Parameter	Value
Product Name	SLC1A1 Rabbit pAb
Remarks/Alias	DCBXA; EAAC1; EAAT3; SCZD18; hEAAC1; SLC1A1
Species	Human
GeneID (Human)	6505
GeneID	6505
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 430-524 of human SLC1A1 (NP_004161.4)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43005
Protein Weight	57kDa
Shipping	Ice bag