

Anti-Slc25A19 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Hc21

Item Number: CM0018113

Introduction

Rabbit polyclonal antibody against human SLC25A19, a mitochondrial thiamine pyrophosphate carrier. Suitable for WB, IF/ICC, ELISA applications with human and mouse reactivity. Useful for studying mitochondrial transport and related disorders.

[Learn More](#)

Parameter	Value
Product Name	SLC25A19 Rabbit pAb
Remarks/Alias	DNC; TPC; MUP1; MCPHA; MTPPT; THMD3; THMD4; SLC25A19
Species	Human
GeneID (Human)	60386
GeneID	60386
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human SLC25A19 (NP_068380.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HC21
Protein Weight	36kDa
Shipping	Ice bag